



Test Certificate

Part **Basket Holder with Stirring Shaft**
Part Number **07558-07**
Serial Number **P13 001085 – P13 001194**



			USP / Ph.Eur. Tolerances
Dimensions	Shaft Ø	9.4 – 10.1 mm	--
	Cap Ø	25.0 mm	± 3.0 mm
	Cap size	5.1 mm	± 0.5 mm
	Clear opening	20.2 mm	± 1.0 mm
	Vent hole Ø	2.0 mm	± 0.5 mm
Material	Stainless Steel (electrolytic polished)		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2023-10-05

MP

05-10-2023 17:07:27

Moreshwar Pulekar
SOTAX Quality Management
C-64 | First Floor | Suchi Business Park | TTC Industrial area | MIDC Turbhe | Navi Mumbai 400705 | India