




# Test Certificate

Part 2 Liter Vessel  
Part Number 08222  
Serial Number 93275 - 93288



		USP / Ph.Eur. Tolerances		
<b>Dimensions</b>	Vessel	Inside Ø a	100 mm	98 – 106 mm
		Height b	200 mm	280 - 300 mm
		With hemispherical bottom		
<b>Material</b>	Borosilicate glass			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2023-03-28

Bruno Balli  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland