




Test Certificate

Part 1 Liter Vessel Inactinic Brown
 Part Number 14143-12
 Serial Number P36 000001, P36 000004, P36 000005, P36 000007,
 P36 000008, P36 000009, P36 000012, P36 000013,
 P36 000014, P36 000015, P36 000017, P36 000024,
 P36 000028, P36 000029.



			USP / Ph.Eur. Tolerances
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm
		Height b	168.2 mm
		With hemispherical bottom	160 – 210 mm
Material	Borosilicate glass inactinic brown		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2021-04-30

4/30/2021 2:59:49 PM

SOTAX Quality Management
 EL 142/143 | Ground Floor | TTC Industrial area | MIDC Mahape | Navi Mumbai 400701 | India