

Test Certificate

Part: **Basket Holder for SOTAX AT MD**
Number: 231647
Serial Number: W1702601- W1702665

For SOTAX AT - Dissolution Testing Systems



Dimensions:	Shaft Ø	9.4 to 10.1 mm	Tolerances
	Cap Ø	25.0 mm	--
	Cap size	5.1 mm	± 3.0 mm
	Vent hole Ø	2.0 mm	± 0.5 mm

Material: 316 Stainless Steel

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. (at time of delivery).



SOTAX Quality Management
Ajay Pydah
June 1, 2018