

# Test Certificate

Part: **Basket, 40-Mesh SS, Dissolution**  
 Number: 63127  
 Serial Number: QLA86591-QLA86649

For SOTAX AT - Dissolution Testing Systems



|             |                       | USP Tolerances |                |
|-------------|-----------------------|----------------|----------------|
| Dimensions: | Basket rim O.D. Ø     | 25.0 mm        | ± 3.0 mm       |
|             | Height                | 37.0 mm        | ± 3.0 mm       |
|             | Wire Ø                | 0.25 mm        | 0.22 – 0.31 mm |
|             | Wire openings         | 0.4 mm         | 0.36 - 0.44 mm |
|             | Screen opening height | 27.0 mm        | ± 1.0 mm       |
|             | Screen O.D. Ø         | 22.2 mm        | ± 1.0 mm       |
|             | Clear opening Ø       | 20.2 mm        | ± 1.0 mm       |

Material: Stainless Steel 316

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management  
 Ajay Pydah  
 July 23, 2018